

CPAP Comfort Cover, LLC

rev. 4.4.17

rev. 4.4.17							
FITS MASK	MASK SIZE	SKUs FOR DIFFERENT COVER SIZES					
		7590	7090	6090	6060	5060	5050
Identifying Color of Dot on Cover		Lime Green	Yellow	Gray	Violet	Light Blue	Turquoise
Fisher & Paykel							
Eson/Eson 2 Nasal **	S-M-L				X		
Flexifit HC431 Full	S	X					
Flexifit HC432 Full	S	X					
Simplus Full **	S		X				
Simplus Full **	M-L	X					
Zest Nasal	Petite-Standard				X		
Zest Q Nasal	Petite-Standard				X		
Lady Zest Q Nasal	Petite-Standard				X		
ResMed							
AirFit F10 Full	XS-S		X				
AirFit F10 Full	M-L	X					
AirFit F20 Full **	S		X				
AirFit F20 Full **	M	X					
AirFit F10/F20 For Her Full	XS-S		X				
AirFit F10/F20 For Her Full	M	X					
AirFit N10/N20 Nasal	Small						X
AirFit N10/N20 Nasal	Standard-Wide					X	
AirFit N10/N20 For Her Nasal	Small						X
Mirage Activa LT Nasal	S-M-L-Large Wide		X				
Mirage Quattro Full *	XS-S-M-L	X					
Mirage Soft Gel Nasal	L-Large Wide		X				
Quattro Air Full	XS-S		X				
Quattro Air Full	M-L	X					
Quattro Air For Her Full	XS-S		X				
Quattro Air For Her Full	M	X					
Quattro FX Full	S-M-L		X				
Quattro FX For Her Full	S-M-L		X				
Respironics							
Amara Gel/Silicone Full	S-M-L		X				
Comfort Gel Blue Full *	S-M-L	X					
Comfort Gel Blue Nasal	Petite				X		
Comfort Gel Blue Nasal	S-M-L			X			
Pico Nasal	S/M-L						X
Pico Nasal	XL				X		
Profile Lite Gel Nasal	L			X			
TrueBlue Nasal	Petite				X		
Wisp Nasal	Petite-Small						X
Wisp Nasal	L-XL				X		

* Loop the cover's elastic over the tabs on the side of the mask to keep it secure.

** Snap the cushion into the mask before attaching a cover.